

NORTHROP EXHIBIT O

M. Ching

Cells of Mito Ching

PCR rxns on miles heaters 35 ² heaters
Volume ~ 25 μ l Volume wells


1) Standard rxns 20, 20, 30, 40 μ l

20, 30, 40 μ l w/ graphite pencil tips

~ 60 μ l oil in each

water temp:
max = 45°C
min = 42°C

2) Device: ~ 60 μ l rxn mixture
~ 200 μ l oil



Cycling w/ 7.7 volts } 1.5 W
200 mA

oil warms?
up

upline = 47 sec / 29 sec = 21st cycle
down = 23 sec / 15 sec = " "

- Thermocouple, touching membrane center
- Type (?)
- Diameter .005" = 100 μ m
tip = 200 μ m

poly = ~~200~~ thickness 0.6 μ m
clutter: 0.3 μ m


Cont

M. Allington

- 30 cycles completed on device

- 25 cycles on standards

gel:





1	20	ul	4 graphite
2	30		
3	40		

4	20	ul	4 out graphite
5	20		
6	20		
7	40		

25
cycles

A - E = device = 30 cycles





no comp

Results:

- Graphite did not have a significant effect (lanes 1-3)
- Primer Dimer formed in wells (lanes 9-11) due probably to not reaching high enough T for lambda to denature
- note this system has 2-base overlap which is least toward primer-dimer formation
- evidence of steep T-gradient
- Try higher T (40°) longer 1 min.