

NORTHROP EXHIBIT H

5/11/92

(1)

Key features:

BSAC

PCR or thermal-cycle driven
synthetic reactions, but
includes other DNA (or
other biochem or chem) reactions

All techns. involving reaction-based
methodologies

Ability to integrate elec, mech
or optical comp's. using micro-fab
techs:

Reaction parameter control: heating,
pumping, circulating, cooling

Reaction or reagent detection-manipulation

Q34 Appl. deifed Fig. 2 chamber &
specific sys. for carrying out
PCR reactions. What are
minimum comp's. of such a system?

Output:
mass change,
density
viscosity

Bio chem reaction manip, control & detection: (1) wave
devices - to pp, mix & detect;
(2) electrokinetic effects to pp,
size of pore; (3) interg-T
control devices; (4) optical filters;
(5) chamber w/ heating elements

ATT

Allen Northrup
Dick White

11/92

→ thermal cycling

claim: mixing chamber; also
clean system.

Micro liter vols → used in system
smaller — can be done faster
w/ less power.

96° ↔ 55°c → heat & cool.
system on

Re-packaged/chips → novel:
integ. on — chips.

* For any other reaction that
requires thermal cycling.
chem. or bio-chem

Monolithic, microFab device →
1st time such a sys
has been done on this scale.

(micro pumps) (3)

use of hand-wheel as γ p. \rightarrow
moves particles in single

* file
* can't see if I have
file or white's u
5,006, 149.

unifex: microheaters in chamber
key part of inv. / system.

Fiber optic inside reaction chamber
for detection.

PCR techs now take about ^{1 to 5 min. not} ~~1 to 2 hr~~
~~at inv. can get results in~~
~~1 to 5 minutes.~~ Fastest

in time. ~ can get
results in a few seconds.

re cycles on the order of
a cycle. Due to small
vols + high surface area.
could be portable + work w/
batteries. as opposed to 110V.

~~to diff. surfaces~~

* Reactions dep. on concentration
at

on \rightarrow 3 volts in 90°C based on heater Res.
 off \rightarrow 1/2 volts in 55°C " "


Monolithic Fab — is very important
 — can do in batches on
 wafers. Thus they are
 cheap & are disposable.

* Need spec. sheet for chamber
 fabrication

Pplng in chamber — re Dick
 white

Chamber — includes ultrasonic
 agitator for mixing.

\downarrow Also useful for cell dispensing



pump, mixer etc run
a diff pres (F1, F2, F3) so
can turn diff cover on top
by opening a diff pres.
Heating is also done by tuning F4
for heater:
Key doing PCR on chip w/
reaction chamber.

~~Materials~~ used in chamber
can negatively effect PCR
reaction eg. presence
divalent cations will prevent
PCR from being.

Reactants on chip in advance
of solid as a unit.

Use of surface tension to maintain reactant


Possible talks / pubs : none on
the horiz; Oct 1992

Draft appl to Investors by beginning
to mid-June & file by end
of June

* Rotate drive

Ads: All parts of channels moved at
same speed - adv. of in-
vent PP.


dee whitened
US 5,006,749



Blast Flow
- more control

Liquid-phase
generators

T^0 control \rightarrow Nice feature.



Use ~~leads~~ ^{electrodes} A + B to measure T^0 change. Send in known current I + measure voltage bet. electrodes A + B.

This detail should be in log.

Si O₂ - is material at bottom
of chamber; sed material
handled by bio-chemists

DNA probes - in chamber;
parts attach to probes + sensor
change so know have attachment
to probes.